





Test kit for manual concentration measurement of pH-neutral cleaning agents

In order to achieve consistent cleaning results, the user has to ensure that the cleaning agent concentration stays within the recommended application range. Factors such as "drag-out", dilution, and evaporation influence the cleaning agent concentration in the cleaning bath, so it is highly recommended to take regular measurements of allprocess related parameters. The ZESTRON® Bath Analyzer 20 is a simple, easy-to-use method which provides reliable and accurate results for either fresh or contaminated cleaning baths.

Simple measurement in only a few steps:



1) Sampling / Decanting



2) Adding test solution



3) Shaking / Evaluating

Advantages of the ZESTRON® Bath Analyzer 20:

- Very precise concentration measurement:
 - » for increasing contamination levels
 - for either or contaminated cleaning baths
- Simple and fast procedure
- Mobile product that can be used anywhere
- Requires little space and minimal investment

Watch the product video on www.zestron.com



Specifically developed for the following cleaning agents:

VIGON® SC 200	VIGON® TC 150	VIGON® PE 180	VIGON® PM 105	HYDRON® SE 230A
VIGON® SC 202	VIGON® N 600	VIGON® PE 200	VIGON® 1000 CR	HYDRON® SC 300
VIGON® SC 210	VIGON® N 640	VIGON® PE 215N	HYDRON® SE 220	ZESTRON® VD 200
ATRON® DC	VIGON® N 680			

The ZESTRON® Bath Analyzer 20 includes:



- 1) PE bottle with injection nozzle
- 2) 2 pair of vinyl gloves
- 3) Digital thermometer (°C/°F)
- 4) Sampling beaker
- 5) Glass cylinder with marking
- 6) 6 x 100 ml ZESTRON® Bath Analyzer 10 test solution
- 7) User manual, including correlation charts (not pictured)

Authorized Distributor Mectronics Marketing Services

313, 2nd Floor, Patparganj Industrial Area,

New Delhi - 110092 (India)

Tel : +91 (011)-43041581 / 42208256

E-Mail: mectronics@mectronics.in

Web: www.mectronics.in